PHM201 Portable pH Meter
PHM210 Standard pH Meter

Operating Instructions
Bedienungsanleitung
Mode d'Emploi

Contents

English

Setting up	3
Starting up	5
Calibration	7
Measurements	9
Printed Messages	11
Troubleshooting	12
Maintenance	14
Connection of Peripherals	15
Remote Control using a PC	16
Accessories	18
Specifications	19

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Introduction

The PHM201 Portable pH Meter and PHM210 Standard pH Meter form part of Radiometer Analytical's MeterLab® range. Ideal for use in the field, the PHM201 packs the features of a benchtop meter into a portable model. Whether used in the lab or for educational purposes, the PHM210 fulfills all requirements for straightforward pH and mV measurements.

Both models automatically recognise buffers used for calibration. Just select which of the three types of buffer you want to use. For automatic temperature measurement, a temperature sensor can be connected. The PHM210 has a built-in interface for connection of printer/PC, and outputs for analog recorder and the SAM7 Sample Stand.

For Operating Instructions in English please turn to page 3.

Turn to the centre pages for your Short-form Instructions; please detach carefully.

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The PHM201 Portable pH Meter

Installing the Battery

Slide open the battery compartment situated on the back panel of the PHM201, in the direction of the arrow. Place a 9 V battery into the compartment and fit the + and - terminals into the corresponding sockets. Close the battery compartment.

Connecting the PHM201

The PHM201 is designed to run on a battery. However, it can be run off the mains by using a mains adapter (optional). Remove the "protective cap" and plug the adapter into the external power input. For convenience, you can keep the "protective cap" in the battery compartment.

Using the mains adapter

230 Vac 50 Hz / 9 V **regulated**, 180 mA 

Attention: we accept no responsibility when using a mains adapter which has not been supplied by Radiometer Analytical.

Always unplug the adapter from the pH meter first and then from the mains socket in order to avoid losing your data.

The PHM210 Standard pH Meter

Connecting the PHM210

Connect the PHM210 to the mains supply via the **POWER IN** (12 V) socket using the AC/DC adapter supplied. For safety reasons make sure that the power adapter is placed at a distance of at least 1 metre from a water supply.

Connecting the Electrode and Temperature Sensor


Mount the electrode on the electrode stand. The E190 Electrode Stand or the SAM7 Sample Stand, which is a combined electrode holder and magnetic stirrer, are ideal for setups with the PHM201/PHM210. The combined electrode is connected to the PHM201 via the BNC socket situated on the right-hand side of the instrument and to the PHM210 via the BNC socket marked **ELECTRODE** on the rear panel of the PHM210. The PHM210 can also use two separate electrodes. Connect the indicator electrode, e.g. pH glass electrode, to the BNC socket and the reference electrode to the Banana socket marked **REF**.

If you are using an electrode mounted with a Type 7 plug with your PHM201/PHM210, connect to the BNC socket via an adapter (part no. A94P802).

If desired, you can use the T201 Temperature Sensor for automatic temperature measurement. This is connected via the CINCH socket on the right-hand side of the PHM201 or to the input marked **TEMP** on the PHM210 rear panel.

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Starting up

1. Switch on your pH meter by pressing the **On** key on the front panel of the PHM201 or the red **POWER ON/OFF** button situated on the rear panel of the PHM210. The first display you see will remain for a few seconds and shows the name of the meter and the software reference while the PHM201/PHM210 is performing an Autotest (see page 4 for an overview of the displays).

Note: hold the **Cal** key down during the Autotest to enter the language display. Using the **left** or **right arrow** key select the language to be used, choose from English, German and French. Press ✓ to continue.

2. The buffer type can now be selected. Use the **left** or **right arrow** key to select **IUPAC Series pH Standards** (pH 1.679, 4.005, 7.000, 10.012 or 12.45), **Technical Series** (pH 1.09, 4.65 or 9.23) or **4-7-10 Series** (pH 4.00, 7.00, 10.00). Press ✓ to continue.

Note: the PHM201/210 display buffer values to 2 decimal places.

3. The PHM201 will now display the **AUTO OFF YES/NO** display. If YES is selected: the PHM201 will switch itself off automatically within 4 minutes if no keys are pressed. Press ✓ to continue.

Note: this display does not appear when using the PHM210.

4. The PHM201/PHM210 now goes to the live pH display. You are now ready to perform a pH calibration or mV measurement.

Calibration Procedure

Display Overview

Last Calibration
pH° 7.00 100.0%✓

Buffers in use 4 seconds
TECHNICAL
IUPAC
4-7-10

Cal Temperature 25.0°C ✓
If no temperature sensor is connected adjust the temperature manually using ◀ ▶

Dip Electrode in Buffer 1 ✓

STA 6.98 pH
Cal -8 mV 21.3°C

Buffer 1: 7.00
Dip in Buf.2 ✓

1 Point Cal
pH° 6.44 100.0%✓

Auto printout of calibration data with PHM210

STA 3.93 pH
Cal 166 mV 21.5°C

STAB 4.01 pH
21.6 °C

With temperature sensor

or

Buffer 2: 4.01
pH° 6.44 99.0%✓

or

STAB 4.01 pH
◀ ▶

Without temperature sensor

Auto printout of calibration data with PHM210

Calibration

A pH meter must be calibrated against standards, i.e. buffer solutions, to adjust the pH reading to the individual characteristics of the electrode used. The PHM201 and the PHM210 can perform either a two-buffer calibration or a one-buffer calibration.

The calibration is performed with automatic recognition of the pH buffer solutions. Three different types of buffer solution can be used: IUPAC Series pH Standards, Technical Series manufactured according to DIN19267 or 4-7-10 Series. Buffers of different types cannot be used in one calibration (e.g. using a Technical buffer followed by a IUPAC standard pH solution for the calibration will lead to erroneous results). The pH value of a buffer changes with temperature; for calibration calculation the PHM201/PHM210 uses the actual pH values at the current temperature. The temperature is entered using the arrow keys in steps of 1°C or measured automatically when a temperature sensor is connected (temperature displayed to 1 decimal place e.g. 23.4°C).

Procedure

Refer to the display overview given on page 6.

1. Switch on the pH meter.
2. Mount and connect the electrodes; refer to page 3.
3. Pour the desired buffer solution into a beaker. Select the buffer close to the pH value of the samples. To obtain the highest accuracy, the pH values of the buffers used during a two-point calibration should lie on either side of the pH value of the samples.
4. Press the **Cal** key to start the calibration. The calibration results obtained during the previous calibration will now appear on the display.
5. Press **✓** and the type of buffer to be used during the calibration is displayed for a few seconds. If you want to use another type of buffer, select this in the setup state; see page 5.
6. If no temperature sensor is connected, adjust the calibration temperature manually using the **left or right arrow** key. Press the **✓** key to start the calibration procedure.
Note: If a temperature sensor is connected, the temperature will be measured directly and the **Cal Temperature** display will not appear.
7. Dip the electrode in the first buffer when the message **Dip Electrode in Buffer 1** appears.
8. Check that the electrode(s) is properly immersed in the buffer solution and that proper stirring is applied. Press **✓** to continue.
9. The calibration's progress can now be followed on the display by means of the pH readings and the visual **STABILITY** indicator.

When the electrode signal is stable the PHM201/PHM210 will show the following display:

Buffer 1: 7.00
Dip in Buf.2 ✓

10. If a one-point calibration is required, press the **pH/mV** key to end the procedure. The following display appears, showing the new zero pH value and the last sensitivity. The PHM210 will automatically print out the calibration data if a printer is connected; refer to page 11. Press ✓ to go to the measurement mode; refer to page 9.

1 Point Cal
pH° 6.92 100.0%✓

11. To continue a two-point calibration, rinse the electrode(s) with deionised water and dip the electrode(s) in buffer 2. Press ✓ to continue.
12. The calibration's progress can now be followed on the display by means of the pH readings and the visual **STABILITY** indicator. When the electrode signal is stable the PHM201/PHM210 will show the calibration result, zero pH and sensitivity. The PHM210 will automatically print out the calibration data if a printer is connected. If an error occurs during the calibration an error message will be printed out; refer to page 11.

Buffer 2: 10.01
pH° 6.92 99.7%✓

13. Press the ✓ key or the **pH/mV** key to go to the measuring mode.

The calibration results are stored in the PHM201/PHM210 memory and used for the pH measurements which follow. If an error message is displayed during the calibration procedure, e.g. **Same buffer change buffer**, refer to "Troubleshooting" page 12.

How to reset the calibration data

The following calibration data are programmed as default values in the PHM201/PHM210: zero pH: 7.00, sensitivity: 100%, temperature: 25 °C. If you wish to return to these values, proceed as follows:

1. Hold down the **Cal** key for 4 seconds.
2. Using the arrow keys select either **Reset Cal Yes** or **No**. Press ✓ to enter.
3. **If Yes:** the **Last calibration** display appears showing the default values. Press ✓ to initiate a new calibration. Go to step 6 of the calibration procedure.
4. **If No:** the **Last calibration** display appears showing the results of the last calibration performed. Press ✓ to initiate a new calibration. Go to step 6 of the calibration procedure.

Measurements


pH measurements can be performed in two ways:

- Manually with a live pH reading on the display.
- Using the AUTOREAD function which locks the result on the display as soon as the electrode signal has stabilised. If a printer is connected to the PHM210 the results are automatically printed out.


Manual pH Measurements

1. Perform a one or two-point calibration as previously described. To check the stored calibration results, press the **Cal** key.
2. Rinse the electrode(s) and temperature sensor thoroughly with deionised water and dab dry with soft tissue.
3. Immerse the electrode(s) and temperature sensor in the sample and apply proper stirring.
4. Press the **pH/mV** key to initiate sample measurements.

If a temperature sensor has been connected, the actual pH value and temperature measured will be displayed. If no temperature sensor is used, the actual measuring temperature should now be adjusted using the arrow keys.

With a temperature sensor

Without a temperature sensor

5. Read the pH value when the signal is stable, i.e. when the stability indicator displays "STAB" continuously.

When "STAB" is displayed, it means that the stability criterion is reached (~ 0.05 pH/min.). The letters " ", "S", "ST", "STA" appear as stabilisation progresses.

pH AUTOREAD

1. Follow steps 1 to 4 on page 9.
2. Press the ✓ key to initiate the AUTOREAD function.
3. The PHM201/PHM210 will now monitor the electrode signal's stability. The measurement can be followed by means of the stability indicator until acceptance occurs.

STA	10.06	pH
AUTOREAD	21.3	°C

4. Once stable, the accepted pH value will be locked on the display and the following display will appear:

RESULT	10.62	pH
	22.1	°C

The PHM210 will automatically print out this result if a printer is connected; refer to page 11.

5. Press the ✓ key to start a new AUTOREAD measurement.

The following display will appear if the ✓ key has been pressed while awaiting stabilisation or if the electrode signal is not stable within 3 minutes. In these situations the result is marked with an asterisk (*).

RESULT	* 11.62	pH
	22.1	°C

The PHM210 will automatically print out the unstable result if a printer is connected; refer to page 11.

Manual mV Measurements

mV measurements are performed in a similar way to pH measurements.

1. Press the **pH/mV** key to display the live pH reading.
2. Press the **pH/mV** key a second time to display the live mV display.
3. Read the mV value when the signal is stable, i.e. when the stability indicator displays "STAB" continuously.

When "STAB" is displayed, it means that the stability criterion is reached (3 mV/min.). The letters " ", "S", "ST", "STA" appear as stabilisation progresses.

4. To return to the live pH display press the **pH/mV** key.

Printed Messages

Calibration:

The following printout is obtained automatically after a two-point calibration:

```
***** PHM210 Standard pH Meter *****
      Calibration
      Temperature  22.1°C
      Buffer 1      10.01
      Buffer 2      7.00
      Zero pH      6.98
      Sensitivity  98.9 %
```

If an error occurs during the calibration procedure, one of the following error messages will be printed out instead of the above calibration result:

```
SENSITIVITY ERROR: 87.1 %
```

or

```
ZERO pH ERROR: 5.75
```

Results:

In pH AUTOREAD measurement mode, the following printout is obtained when the accepted pH value is locked on the display:

```
No. 1      7.56 pH  21.7°C
No. 2      7.56 pH  22°C
No. 3      *6.93 pH 21.7°C
up to...
No. 9999  *10.34pH 23°C
```

Note: If the pH value is unstable the result is marked with an asterisk (*). The result number returns to No. 1 each time a new calibration is performed, when the PHM210 is switched off or when No. 9999 has been reached.

Troubleshooting

A number of tests are automatically performed when the PHM201/PHM210 is switched on, during calibration and during measurement. If an error is detected, an error message in clear text will appear on the display. The messages are listed below.

Error Messages

Message	Error/Action
BATTERY LOW	The stored data are still reliable. Change the battery, refer to page 14.
REPLACE BATTERY USER DATA LOST	The stored data have been replaced by default values. Change the battery, refer to page 14.
USER DATA LOST	The stored data have been replaced by default values when the battery in the PHM201 has been removed or the adapter has been unplugged from the mains socket; see page 3.
OUT OF RANGE	This message appears if the measured potential exceeds the -1999 to +1999 mV or -9 to +23 pH range. Check electrode condition and electrode connection.
UNKNOWN BUFFER REPEAT MEAS.	The buffer is not recognised. Check that the buffer used is of the selected type, and that the solution is fresh. Select your buffer type in the setup state; see page 5. Check electrode condition and electrode connection before repeating the measurement.
UNSTABLE SIGNAL REPEAT MEAS.	The electrode signal has not stabilised within 3 minutes during calibration or during an AUTOREAD measurement. Check electrode connections, stirring conditions, etc. Press the ✓ key to repeat measurement or the pH/mV key to return to manual measurement mode.

Message	Error/Action
SAME BUFFER CHANGE BUFFER	<p>The same buffer has been used twice for a two-point calibration.</p> <p>Calibration against the correct second buffer can be initiated after pressing the ✓ key.</p>
SAME BUFFER OR ELECTRODE ERROR	<p>This message will appear during a two-point calibration if a pH 7 buffer solution was used as buffer 1 and buffer 2 or if the electrode is not functioning properly.</p> <p>Check the electrode condition and make sure the buffer has not been used twice. Press the ✓ key to restart the calibration in buffer 2.</p>
ZERO pH ERROR	<p>This message appears if the zero pH is less than 5.8 or greater than 7.5. If the PHM210 is connected to a printer this message is automatically printed out.</p> <p>Check that the correct buffer has been used. Check the electrode.</p>
SENSITIVITY ERROR	<p>The electrode sensitivity has not been accepted as it is outside the 95 to 102 % range. If the PHM210 is connected to a printer this message is automatically printed out.</p> <p>Check the electrode connection and condition.</p> <p>Press the ✓ key to restart the calibration in buffer 2.</p>
PRINTER ERROR	<p>If a printer is connected to the PHM210 and data transmission cannot be performed properly, this message will appear.</p> <p>Check that the printer has been turned on and that the paper is in place.</p>
CALL FOR SERVICE Test 1 failure	<p>The battery has not been replaced immediately. Change the battery, refer to page 14.</p> <p>If after replacing the battery the error persists, contact your Radiometer Analytical Representative.</p>

Note:

In the unlikely event of electrostatic energy being discharged onto the PHM210, interference may occur to the display. This is a temporary situation - the display will return to normal within a few seconds.

Maintenance

The PHM201 and the PHM210 require minimum maintenance. However, the outer surface of the pH meters should be regularly cleaned with a soft cloth and tepid water. No strong acids or bases nor organic solvents should be used to clean the pH meter.

Changing the Battery in the PHM201

The PHM201 runs on a 9 V battery. The battery should be replaced when the warning message telling the user that the battery is running low appears; see page 12.

It is recommended to change the battery as soon as possible to avoid the loss of calibration data.

To change the battery, slide open the battery compartment, situated on the back panel of the PHM201 in the direction of the arrow, take out the old battery and replace it with one of the same type and rating (9 V). Close the battery compartment.

Adjusting the Contrast

The contrast of the PHM201's display can be adjusted by means of a screw in the battery compartment. By turning the screw clockwise or anticlockwise, the display will turn darker or brighter respectively.

RECOMMENDATION:

Unplug the electrodes before switching off your pH meter.

Connection of Peripherals


The PHM210 is equipped with terminals for connection of Printer/Personal computer, recorder and the SAM7 Sample Stand. The terminals on the rear panel are illustrated on page 2.

PRINTER/PC:

For the connection of a printer or a personal computer.

RS232C specifications:

9-pin, sub D-connector, (9600, n, 8, 1): 9600 baud, noparity, 8 data bits, 1 stop bit.


Pin connections:

- | | |
|----------------------------------|-----|
| 1. No connection | |
| 2. Receiving Data (input): | RxD |
| 3. Transmitted Data (output): | TxD |
| 4. Data Terminal Ready (output): | DTR |
| 5. Ground: | GND |
| 6. No connection: | DSR |
| 7. Request to Send (output): | RTS |
| 8. Clear to Send (input): | CTS |
| 9. No connection | |

RECORDER:

8-pin DIN

The PHM210 offers two types of analogue output signals:


Pin connections:

1. Direct electrode signal output
2. Calibrated pH output
1 pH unit = 100 mV
Resolution = 1 mV
5. 0 V

POWER OUT:

6-pin, DIN

The PHM210 can supply power to the SAM7.

Pin connections:

3. Protective ground
6. -12 V, 0.3 A output

Remote Control using a PC

Communication protocol:

The PHM210 can be controlled from a PC.

Data transmission: (9600, e, 8, 1): 9600 baud, no parity, 8 data bits, 1 stop bit.

Description of messages:

Command	Answer	Comments
?		List of commands
?TY	=TY PHM210 X9543	PHM210 software reference
?ST	=ST mV =ST pH =ST CAL =ST TEST mode	Status: mV measurement Status: pH measurement Status: Calibration mode Status: Autotest and setup
?RD t	=RD xx (.x) C	Read temperature
?RD mV	=RD xxxx mV	Read potential
?RD pH	=RD xx.xx pH	Read pH
?MR adr	=MR val	Memory Read address value 1 E ₀ (in mV) 2 S ₂₅ (in mV) 3 Type of buffers 0: (TECHNICAL) 1: (IUPAC) 2: (4-7-10) 4 Language E: (ENGLISH) D: (DEUTSCH) F: (FRANCAIS)

Command	Answer	Comments
>MW adr, val	<MW	Memory Write address value 1 E' ₀ (in mV) 2 S ₂₅ (in mV) 3 Type of buffers 0: (TECHNICAL) 1: (IUPAC) 2: (4-7-10) 4 Language E: (ENGLISH) D: (DEUTSCH) F: (FRANCAIS)
>MU	<MU	Updates the contents of the non volatile memory

Keyboard emulation:

Command	Answer	Comments
}S	{S	pH/mV key
}C	{C	Cal key
}O	{O	✓ key
}R	{R	► key
}L	{L	◄ key

Accessories

The PHM201 Portable pH Meter (part no. R21M001) is supplied in a handy slimline carrying case containing three measuring beakers and three empty bottles for buffers with room for an electrode (not supplied).

One of the following electrodes is suggested for use with the PHM201/PHM210:

pHC3001-8 Combined pH Electrode	E16M300
pHC3005-8 Combined pH Electrode in robust epoxy body	E16M302
pHC3105-8 gel-filled Combined pH Electrode in robust epoxy body	E16M304
pHC3081-8 Combined pH Electrode with built-in temperature sensor	E16M305
pHC2085-8 Combined pH Electrode in robust epoxy body with a built-in temperature sensor	E16M501

Accessories Available

T201 Temperature Sensor	E51M001
Solution Kit: Buffer solutions: pH 4, pH 7 and pH 10, 500 ml of each. 100 ml KCl/AgCl solution.	S91M003
Power adapter for the PHM201 230 Vac 50 Hz / Regulated 9 Vdc, 180 mA	A66B003
9 V Battery (PHM201)	A65B001
Serial cable 9-25 pin (PHM210 to PC)	A95X502
Serial cable 9-9 pin (PHM210 to PC)	A95X501
Serial cable 9-25 pin (PHM210 to printer)	A95P201
Recorder cable	A95R102
E190 Electrode Stand	809-206
SAM7 Sample Stand 230 V version	391-563
115 V version	391-564
Connection cable (PHM210 to SAM7)	617-896

Specifications

Measuring ranges:	pH : -9.00 to 23.00 mV : -1999 to +1999 °C : 0 to +99.9
Resolution:	pH : 0.01 mV : 1 on the display °C : 1 with no temp. sensor 0.1 with temp. sensor
Accuracy:	pH : ±0.02 pH after calibration mV : typical value: ±0.1% of reading or ±1 mV (whichever is greater)* °C : ±0.5°C
Electrode Inputs:	BNC plug for combined or glass electrode CINCH plug for temperature sensor PHM210: Banana plug for reference electrode
Measurement procedures:	Manual pH and mV reading. AUTOREAD pH: result is locked on display when stability criterion is reached within 3 min.
Stability criterion:	3 mV/min. ~0.05 pH/min.
Calibration:	One or two-point calibration
Auto recognition of buffers:	IUPAC Series pH Standards: pH 1.679, 4.005, 7.000, 10.012 and 12.45 Technical Series (DIN19267): pH 1.09, 4.65 and 9.23 4-7-10 Series: pH 4.00, 7.00 and 10.00
Electrode requirements:	Sensitivity: 95 to 102% Zero pH: 5.80 to 7.50
Electrode input resistance:	> 10 ¹² Ω
PHM210 terminals:	RS232C port for printer or PC MAB6 - for supplying power to the SAM7 MAB8 - for potential and pH recorder outputs
Finish:	Splashproof casing 2 x 16-character alphanumeric LCD display.
Languages:	English, German and French
Ambient temperature:	5 to 40°C
Relative humidity:	20 to 80%
Power requirements - PHM201:	9 V alkaline battery or mains adapter
Power requirements - PHM210:	12 V mains adapter
Weight - PHM201:	0.4 kg
Weight - PHM210:	1.5 kg

*mV : max. value: ± (0.1% FSR + 1 resolution digit) mV

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Printed by Radiometer Analytical S.A. • France • 2001-10G