
Labofuge[®] 200

SERVICE MANUAL

P/N 12003630

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TABLE OF CONTENTS

Sect.	Title	Doc.- No.	Page
1	OPERATING INSTRUCTIONS (not part of this manual)		
1.1	Brochure	Labo200_E	
2	SERVICE	Labo200_2	
2.1	Servicing Schedule	" - "	2-1/2
2.2	Trouble Shooting	" - "	2-3/5
2.3	Test Points	" - "	2-6
2.4	Cleaning of Instrument Parts	" - "	2-7
2.5	Electrical Safety Check	" - "	2-8
3	FUNCTIONAL DESCRIPTION	Labo200_3	
3.1	Block Functions	" - "	3-1
3.2	Functions of Main Board	" - "	3-1/2
4	DIAGRAMS	Labo200_4	
4.1	Block Diagram	" - "	4-1
4.2	Wiring Diagrams	" - "	4-2/3
4.3	Wiring Connection Diagrams	" - "	4-4
4.4	Main Board 28 230V - Component Plan	" - "	4-5
4.5	Main Board 28 230V - Wiring Diagram	" - "	4-6/8
4.6	Main Board 178 120V - Component Plan	" - "	4-9
4.7	Main Board 178 120V - Wiring Diagram	" - "	4-10/11
5	DISASSEMBLY OF INSTRUMENT PARTS	Labo200_5	
5.1	Housing / Casing Parts	" - "	5-1
5.2	Electrical Components	" - "	5-2
5.3	Drive Components	" - "	5-3
5.4	Mechanical Components	" - "	5-4
5.5	Service Kit	" - "	5-5
6	SPARE PART FIGURES AND LISTS		
6.1	Break Down Drawings	Labo200_6	6-1/2
6.2	Spare Part Lists		
7	Preventive maintenance checklist		
8	TECHNICAL BULLETINS		

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LABOFUGE® 200 COMPACT CENTRIFUGE

Easy to Use

The Heraeus® Labofuge 200 is functional and sophisticated with an ergonomic design. It is ideal for use in medical practices, clinical and small laboratories and as a stand-by unit in large laboratories.

Benefits

- Easy to Use
- Able to run wide range of most popular medical tubes
- Maintenance free
- Safe

TECHNICAL DATA

Set the required speed...

...and set the required time

Labofuge 200	230 V version	110-120 V version
Description	Microprocessor controlled table top centrifuge	
Maximum speed (rpm)	5300	5000
Minimum speed (rpm)	1600	1600
Maximum RCF (x g)	3030	2600
Minimum RCF (x g)	270	270
Maximum capacity (ml)	12 x 15 (glass)	12 x 15 (glass)
Controller	Microprocessor controller for time and speed	
Drive	Brushless induction motor, microprocessor controlled	
Runtime	1-99 min. and continuous operation (Hold)	
Program memory	Stores values last entered, unlimited	
Safety	Lid lock and lid interlock, steel armoured guard ring	
Design	Fibre glass reinforced polyamide housing with high noise and vibration insulation properties	
Dimensions (H x W x D) (mm / inch)	240 x 284 x 375 / 9.5 x 11.2 x 14.8	
Weight (incl. rotor) (kg / lbs)	Approx. 10.7 / 23.5	Approx. 10.7 / 23.5
Power consumption (W)	65	65

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2.1 Servicing Schedule (yearly procedure recommended)

2.1.1 Maintenance Routine without Dismantling the Centrifuge

2.1.1.1 Electrical Installation and Safety

- Switch OFF the centrifuge and disconnect the unit from power, check voltage supply and mains fusing (16 Amps, slow blow characteristic)
- Check condition of plug and wall socket - (let) replace defective parts
- Check cord condition and fixing / connection - replace or refit it

2.1.1.2 Location and Mechanical Installation

- Check the base (ground, table, lorry with lockable wheels etc.) For resonance-free and stable conditions
- Check for a well ventilated place and sufficient distances to walls or adjacent equipment, without exposure to direct sunlight
- Check the leveling of instrument (use a spirit level)

2.1.1.3 Lid Locking Mechanism and Safety

- Connect the centrifuge to power and switch ON
- Check for easy lid closing and self-acting lid opening - if in disorder, readjust lid's swivel hinge and/or lock assembly
- Check the central rubber gasket for lid sealing and replace it, if damaged
- For checking the safety: start the centrifuge, let it shortly run and stop it, the lid must not be unlocked by the microprocessor until the "end" message is shown on the display - if safety circuit is out of function, replace the main board

2.1.1.4 Cleanliness of Spin Chamber and Motor Casing

- Open the lid and remove the rotor (for loosening turn the central winged nut anti clockwise from the motor shaft)
- Clean the spin chamber with a dry and absorbent cloth (remove all dust and moisture - see also 2.5 Cleaning)
- Check the cleanliness of the motor casing and take care of the annular slot around the motor shaft: **penetrating fluids can damage the upper motor bearing**, remove fluids with an injector and/or absorbent paper

2.1.1.5 Rotor and Accessories Condition and Sealing

- Check the condition of rotors and accessory parts (especially all supporting or stressed partitions): the rotor and/or accessory parts must not be used any longer, if there are visible traces of mechanical damage

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2.1.1.6 Rotor Fixing and Motor Shaft

- Check the trouble-free condition of the locking nut and replace it in case of malfunction
- High performance rotors made of plastic have a limited service life. Due to safety reasons they have to be exchanged after 10,000 cycles or after being in use for five years. The expiry date can be seen on the date clock in the rotor lower shell. The two-digit number in the middle represents the year and the arrow points to the expiry month. The number of cycles may be checked at any time by closing the lid and keeping both “set” keys pressed. The accumulated rotor cycles will be displayed in the speed display.
- Check the condition of the drive motor shaft: the centrifuge must not be used any longer, if the drive shaft is damaged (bend, thread is worn out, horizontal grooves etc.)

2.1.1.7 Temperature Level

- Check the air inlet underneath the lid and under the bottom plate for free ventilation, insufficient air flow will lead to temperature rise of rotor, motor and electronic parts

2.1.2 Maintenance Routine after Dismantling the Centrifuge Casing**2.1.2.1 Motor Supporting Elements**

- Check the supporting and damping elements of the drive motor and replace them in case of increased rubber abrasion or abundance of imbalance but at least every 3 years

2.1.2.2 Braking Circuit

- Check the function of the braking circuit (even brake effect)

2.1.2.3 Lead and Screwing Connection

- Check the terminal and plug connections of all leads and on all boards and electrical components, tighten all loosen screwing connections, refit or replace defective parts
- Check the screwing connections of all boards, mechanical and electrical components and re-tighten them if necessary (use screw locking lac for motor mounts and lid lock assembly)

2.1.2.4 Protection Earth Core and Grounding Connections

- Check the protection earth core for continuity and all grounding plug connectors
- Check isolation resistance and accessible current (see 2.5)

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2.2 Trouble Shooting

Error Indication	Error Cause	Possible Error Source	Corrective Procedure
Displays remain dark	No mains voltage supply	Mains fuse or circuit breaker failed	Check fuse or circuit breaker, replace or switch on again
		Defective mains cord	Check instrument cord, replace defective parts
		Defective unit fuse or fuses on main board	Replace it, if fuse blows again, disconnect electrical parts, otherwise replace main board
	No low voltage supply	Faulty connection from indication to power board	Check connections on CPU, indicat. Board and connecting leads, replace defective parts
		Faulty indication or power board	Replace main board
Displays are illuminated, but drive doesn't start	Motor over-temperature switch has tripped	Motor temperature is higher than 120°C	Let motor cool down, then check temperature switch and leads with Ohmmeter
	Rotor didn't turn	Rotor is jammed	Check for easy rotor movement, remove any jamming objects
		Motor is jammed	Replace motor
	Motor didn't start	Connections inter drive and main board	Check terminal and lead connections, replace faulty parts
		Defective drive	Check resistance of motor windings, replace faulty parts
		Faulty condenser	Replace condenser
		Faulty main board	Remove main board completely and replace it
Drive makes noises-no good separation result	Mechanics	Wear out of motor rubber mount	Replace motor rubber mounts
		Motor bearing	Change motor completely
	Electrical	Defective terminal connecton, faulty lead or motor winding	Check voltage on motor terminal and winding resistances -see test points on boards
		Defective driving	Replace main board

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Trouble Shooting

Error Indication	Error Cause	Possible Error Source	Corrective Procedure
Drive doesn't decelerate	No brake current	Faulty main board	Remove main board completely and replace it
Lid cannot be opened by key pressure at standstill	Lid coil is not or not sufficiently supplied with voltage	Missing mains voltage	Remedy see above, manual opening only at standstill
		PTC resistor has released	After a waiting time of 1-2 minutes press key again
		Faulty driving or triac circuit	Replace the complete main board
	Faulty lid coil	Faulty winding of coil	Replace complete lid lock
Lid is not correctly locked	Lid bolt is jamming	Push lid into lock and press the key again	
	Lid is deformed or disadjusted	Readjust the lid centrally	
„Lid“ message appears in speed display	Lid was opened manually during run	Forbidden intervention Emergency opening can only be used at standstill	Close lid immediately, turn power off/on, wait for termination of br phase until end message appears
	Protection circuit (15V) for lid control was interrupted during the run	Defective micro switch or leads or connectors to micro switch are interrupted	Check leads and connectors to micro switch, in case of a faulty micro switch, replace lid lock device completely
„OPEN“ message in speed display	15V supply circuit is interrupted at standstill	Defective micro switch or leads or connectors to micro switch are interrupted	Check leads and micro switch
„br“ message appears in speed display	Rotor comes to standstill without braking force	Short interruption of mains supply	Wait for rotor standstill (appr. 75 seconds) and re-start

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Trouble Shooting

Error Indication	Error Cause	Possible Error Source	Corrective Procedure
„E-“, „2“ message appears in speed field	Maintenance counter expired	Rotor has reached 10000 cycles	Replace rotor, NV-Ram and motor supports (Service Kit)
„E-“, „3“ message appears in speed field	No brake current	Faulty main board	Remove main board completely and replace it
„E-“, „12“ message appears in speed field	Checksum error of NV- RAM	NV-RAM is not initialized or false	Check NV-RAM and socket, insert the correct NV-RAM
	Disturbed data transfer from NV- RAM	Faulty main board	Remove main board completely and replace it
All dots in the display are illuminated	Maintenance counter expired	Rotor has reached 10000 cycles	Replace rotor, NV-Ram and motor supports (Service Kit)

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2.3 Test Points

Test Points	Unit value	Conditions
Mains terminal XA	230V AC 120V AV	All given voltage and current values refer either to 230V ($\pm 10\%$) or 120V ($\pm 10\%$)
Lid micro switch Plug XB	320V DC 160V DC	Voltage drop by open lid at 230V units Voltage drop by open lid at 120V units
Terminal XC Motor voltage	approx. 58V AC 105V AC 120V AC 80V AC 75V AC	Rotor #3760 not loaded, in each case measured after reaching the selected speed 1600rpm 4000rpm 5300rpm at the beginning of the braking phase at the end of the braking phase
Motor current I_M	approx. 1,2A 1,2A 0,8A 3,3A 2,0A	Rotor #3760 not loaded, in each case measured after reaching the selected speed with a soft iron or digital effective measuring instrument 1600rpm 4000rpm 5300rpm maximum during acceleration maximum during braking phase
Motor windings resistance 20°C -insulation value	14,6Ω 14,6Ω > 10MΩ	switch OFF unit, pull off motor plug mains winding (black – green) auxiliary winding (black – yellow) resistance inter each phase and motor casing
Lid solenoid terminal XD	105Ω 29Ω	switch OFF unit, resistance at 20°C (68°F) at 230V units resistance at 20°C (68°F) at 120V units

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2.4 Cleaning of Instrument Parts

ATTENTION - WARNING!

The electrical and electronic components must not be cleaned with moist detergents!

For Cleaning the centrifuge housing or its accessories see Operating Instructions section 1 (maintenance and care).

- **Electronic components**

Clean dusty components carefully with a dry and soft brush and remove loose dust with a vacuum cleaner

- **Vent holes**

Remove dirt from vent holes of front panel or bottom plate using a brush and vacuum cleaner

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2.5 Electrical Safety Check

ATTENTION!

A final electrical safety check must be performed after each maintenance and/or repair!

- **Resistance check of protective conductor**


The measuring value of the resistance between the mains plug's grounding pin and the grounding conductors of the motor, electronic chassis and the casing must not exceed 200 mΩ.

- **Insulation resistance Check**

Check also the insulation resistance between the poles of the mains plug and the grounding conductor; the resistance value must be more than 2 MΩ.

- **Accessible current measured to EN 61 010**

The accessible current must not exceed 3.5 mAmps in single fault condition (interrupted protection earth wire)! In accordance with the EN61010, IEC1010 and UL3101 such a fault condition can be reproduced by the following measuring circuit.

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3.1 Block Functions

The Labofuge 200 is a microprocessor controlled laboratory tabletop centrifuge with induction drive motor and integrated air cooling system.

The unit incorporates following boards and components (see block diagram 4-1):

- Main board 28 (for 230 V mains supply) or main board 178 (for 120 V mains supply) with microprocessor part and power electronics
- Key and indication board (programming: MEGACONTROL simple), part of the main board
- 1 phase induction motor with phase shift capacitor and integrated thermal over temperature switch (C. O. 125°C)
- Double lid lock assembly with solenoid and integrated micro switch (mechanical bolt keeper, magnetically dislocking), attached in front on both sides of the vessel

3.2 Main Board Functions

The main board 28 / 178 is mounted completely behind the front panel. The components on main board are arranged in following groups (see wiring diagram page 4-2).

- Board fusing (2 x 6.25 / 6,3 Amps. slow blow feature)
- Noise filter in accordance with EN 55011
- Power pack for low voltage supply of microprocessor part and power electronics (potential separated by transformer)
- Triac control circuit for lid solenoid (DC supply via bridge rectifier)
- Bridge rectifier for DC intermediate circuit supplying brake path and frequency converter
- Microprocessor part with CPU
- Exchangeable NV-RAM containing specific data of the unit

3.2.1 Power Pack

The power pack (transformer, bridge rectifier and voltage regulator) generates:
U1 = 5V: supplies central processor, key and indication board
reference potential: protective conductor (GND)

3.2.2 Intermediate Circuit with Brake Path and Frequency Converter

The DC intermediate circuit serves as an energy store between the AC power input and the transmitted motor performance. The intermediate circuit consists of a bridge rectifier (4 diodes) and two serial connected reservoir capacitors. The intermediate circuit is over voltage protected.

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Brake Path

To braking the centrifugal drive the induction motor is supplied with DC pulses from the intermediate circuit (affects as frequency zero). The pulse frequency is variable in dependence of the braking time (equal to speed), as a result of which the noise level is changed.

Frequency Converter


The frequency converter provides the drive motor with squarewaved voltage blocks having an amplitude of the half of the DC intermediate circuit voltage. These blocks are variable controlled in frequency and pulse-width modulation. Beginning with low frequency and small pulse-width blocks, both parameters will be increased up to the operators set speed. The set speed value will be reached with a tolerance of $\pm 5\%$ (within the valid mains voltage range)


3.2.3 Micro controller (Central Processing Unit)


The controller block includes the central processor unit (CPU) 87C51 and the data storage (NV-RAM) 9346. The ROM which contains the controlling program is integrated in the CPU. The most important operating parameters (e.g. the maximum speed or the last operator settings) are stored in the non volatile (NV) -RAM.


3.2.4 Key and Indication Board


The key and indication board is part of the main board. The connection to the main board is done by a 14 polar flat cable. The 7-segment displays, the control key LED's and the operating keys are managed by the CPU in multiplexing process. For the setting procedures serves MEGACONTROL "simple" (see Instruction Manual)

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





The index numbers stated in () reappear within the breakdown drawings and the spare part lists.

5.1 Dismantling the Housing

5.1.1 Front panel (102)

- Pull out the mains plug and place the instrument to the rear for the wide of the front panel
- Unscrew both screws (type Ejot 101) at the lower rim of front panel (use a small Phillips screw driver)
- Pull the panel's lower rim a little to front (but no more than 15° to avoid the braking off of the plastic hooks at the upper rim)
- Push the front panel in this position to the top to unlock it's casing connection and depose the panel in front of the unit
- Reassemble the front panel analogously in reverse order and press the lower rim tight to the casing when refitting the Ejot screws (but don't tighten the Ejot screws too much to avoid an overscrewing of the casing's plastic thread)

5.1.2 Lid (210)

- Open the lid and pull out the mains plug
- Unscrew both threaded hinge bolts (214), remove them together with bushings (212), spring washers (213) and rubber O-rings (215) and take off the lid
- Reassemble the lid analogously in reverse order (don't tighten the bolts too much, to avoid a squeezing off of the O-rings)

5.1.3 Dismantling the Casing (205)

- Open the lid and pull out the mains plug
- Unscrew the rotor and remove it from the motor shaft
- Dismantle front panel and lid - see 5.1.1 and 5.1.2
- Remove both Phillips screws (204) from bottom plate in front and lateral from armored chamber
- Disconnect the mains cable from terminal **XA** on main board and pull of the grounding plug
- Pick up the front panel, turn it to 90° and place it in front of armored chamber
- Press the casing upwards and remove it from the bottom plate passing the front panel in vertical position
- Reassemble the instrument analogously in reverse order

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5.2 Replacement of Electrical Components

5.2.1 Fuses on Main Board (113) replacement

- Dismantle the front panel (see 5.1.1)
- Check the micro fuses and if necessary replace them by equivalent and undamaged parts
- Perform a test run and search for blowing cause
- Reassemble the device analogously in reverse order

5.2.2 NV-RAM on Main Board (112)

- Dismantle the front panel (see 5.1.1)
- ATTENTION - CMOS components! Discharge your body before handling! Notice correct position of NV-RAM (8 pins) and pull carefully out of socket
- Reinsert the new NV-RAM correctly
- Reassemble the device analogously in reverse order reconnect it to power and perform a test run

5.2.3 Main Board Replacement (110)

ATTENTION! The key and indication board is part of the main board and cannot be ordered separately!

- Dismantle the front panel (see 5.1.1)
- Disconnect plug connections for lid switches **XB**, protecting earth conductor **PE**, and unscrew the lines for lid solenoid **XD**, motor **XC** and mains supply **XA**
- Remove all screwing of the main and indication board (Ejot type 111 and 122) with fiber washers (121) and take out both boards
- Touch a grounded receptacle to discharge your body before touching the sensitive CMOS components! Take a new NV-RAM (112) out of box (or re-use the old but trouble-free component) and insert it into the socket of the new main board (if necessary, remove a placed but non programmed NV-RAM before)
- Place the parts of the new board correctly and mount it with all screws and fibre washers
- Before assembling the front panel check the trigger function of all keys and correct the board mounting if necessary)
- Do not mix up disconnected cables during re-connection!
- Reassemble the device in reverse order and perform a test run, making sure the drive turns in the right direction (see imprinted arrow on rim of rotor chamber)!

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5.3 Replacement of Drive Components

5.3.1 Disassembly of Drive Motor (310)

- Remove the casing (see 5.1.3)
- Disconnect the motor leads from terminal **XC** and pull off the plug connectors from the capacitor (327)
- Take the leads out of the cable holder (308)
- Unscrew the 3 nuts (313) and remove them together with the lock washers (314)
- Lift the motor out of the rubber mounts and remove the motor from the bottom plate
- Reassemble the motor analogously in reverse order and tighten the 3 nuts with use of torque key (5Nm) and secure them with locking lac
- Reassemble the device in reverse order
- Perform a test run and check the sense of rotation (see imprinted arrow direction on casing)

5.3.2 Motor supports (316)

All 3 rubber supports (316) have to be replaced at the same time and at least every three years! By the way the 6 double-sided sandpaper discs (315) must be exchanged, too.

- Unscrew the drive motor (see 5.3.1), hold it with one hand and lay the armored chamber together with motor onto the side
- Pull the motor from the supports and lay it carefully onto the side
- Remove the 3 screws (224) from the bottom plate and take out the rubber supports together with sand paper discs (315)
- Install the new rubber mounts together with new sandpaper disks, put the lock washers in place, tighten the screws evenly and secure them with screw locking lac
- Put the unit again onto it's feet, mount sandpaper disc, motor and lock washer and tighten the 3 nuts evenly by use of a torque key (5 Nm) and secure them with screw locking lac
- Reassemble the device in reverse order
- At last perform a test run

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5.4 Replacement of Mechanical Components

5.4.1 Lid Lock Assemblies (300)

ATTENTION!

Both lid lock assemblies inclusive micro switches are combined to one spare part and cannot be ordered separately!

- Dismantle the casing (see 5.1.3)
- Disconnect the leads of the lid solenoids from terminal **XD** and pull out the plug **XB** for the micro switches
- Unscrew both upper screws (305) of the lid lock assemblies' attachment from the inside of the armored chamber and the lower Ejoy screws (304) from the outer side
- Push down the levers of micro switches, guide the micro switches through the openings of the armored chamber and remove both lid lock assemblies completely with rubber base plates
- Reinstall the new lid lock assemblies analogously in reverse order and secure the thread of all upper fixing screws (305) with Loctite 221!
- Reassemble the device in reverse order
- Check the lid's correct closing and self-acting opening functions after the correct installation and replace the rubber profile (206)

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
5.5 Service Kit for Labofuge 200 / Medifuge 200


This service kit includes all parts which have to be replaced after termination of 10000 cycles. The kit includes different NV-RAMs which are destined for the use of the referring control program (signification of micro controller No.) on the specific main board.

ATTENTION! When changing the NV-RAM it is absolutely necessary to insert the right one into the main board with the right micro controller No.! Misuse can lead to main board destruction!

The specific variants can be identified as follows:

Cat. No. of unit	Mains connection	NV-RAM signification	Controller signification	Additional identifying information and further notes
3630 3476 3627 15258	220V/50Hz 220V/50Hz 220V/50Hz 220V/60Hz	Labo 8 or Labo ALT	0436	minimum speed setting 3000 min-1, E-02 message, main board #150043 (components on both sides)
3630 3476	220V/50Hz 220V/50Hz	1052 V01	0465	minimum speed setting 1600 min-1, no E-02 message, main board #150043 (components on both sides) or main board #150061 (SMD)
3630	220V/50Hz	1052 V03	0474	cycle counter is set to 10500 the first time, main board 28 #150103 (SMD) or main board #150061 (older version V02 is cancelled)
3631	120V/50Hz	1066 V02	0474	main board 14 #150076 or main board 178 #150178

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Kendro service Information

Heraeus Centrifuges

PREVENTIVE MAINTENANCE CHECKLIST

Labofuge 200

Medifuge 200

Type

S/N

Rotor Inspection

- o Visual inspection
- o Rotor locking screw
- o Expiry date
.....
- o Cycle counter checkup
.....

Preliminary Checks

- o Cabinet
(clean, undamaged)
- o Rotor chamber
(clean)
- o Key foil
(clean, undamaged)
- o Installation environment
(level & clearance)

Performance Checks

- o Acceleration time _____seconds
- o Instrument top speed _____rpm
- o Deceleration time _____seconds
- o Time (10 minute set) _____actual min:sec

Safety test

- o Accessible leakage current <3500 μ A
- o Insulation resistance >10 Megaohm
- o Earth Conductor Resistance <0.2 Ohm

Kendro service Information

Comments: